

Anti-Ano1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Xxa6

Item Number: CM0014823

Introduction

Rabbit polyclonal ANO1 antibody validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Target is a 114kDa calcium-activated chloride channel (TMEM16A/DOG1) critical for ion transport and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	ANO1 Rabbit pAb
Remarks/Alias	DOG1, INDMS, TAOS2, ORAOV2, TMEM16A, ANO1
Species	Human
Gene ID (Human)	55107
Gene ID	55107
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of human ANO1 (NP_060513.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5XXA6
Protein Weight	114kDa
Shipping	Ice bag