

Anti-Ano2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nq90

Item Number: CM0029993

Introduction

Rabbit polyclonal antibody targeting human ANO2 (Anoctamin-2/TMEM16B), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Involved in calcium-activated chloride channel activity in olfaction and vision.

[Learn More](#)

Parameter	Value
Product Name	ANO2 Rabbit pAb
Remarks/Alias	C12orf3; TMEM16B; ANO2
Species	Human
Gene ID (Human)	57101
Gene ID	57101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ANO2 (NP_001265526.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ90
Protein Weight	114kDa
Shipping	Ice bag