

Anti-Anp32A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P39687

Item Number: CM0028256



Introduction

ANP32A Rabbit pAb validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Involved in tumor suppression, apoptosis, cell cycle, transcription. Research use in chromatin and virology.

[Learn More](#)

Parameter	Value
Product Name	ANP32A Rabbit pAb
Remarks/Alias	LANP; MAPM; PP32; HPPCn; PHAP1; PHAPI; I1PP2A; C15orf1; ANP32A
Species	Human
Gene ID (Human)	8125
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39687
Protein Weight	29kDa
Shipping	Ice bag