


Anti-Anp32B Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q92688

Item Number: CM0005458



Introduction

Designed for WB, IF, IHC, ELISA with human, mouse, rat ANP32B. Rabbit monoclonal antibody targeting the multifunctional protein involved in apoptosis, cell cycle, transcription, and viral replication.

[Learn More](#)

Parameter	Specification
Product Name	ANP32B Rabbit mAb
Remarks/Alias	APRIL; SSP29; PHAPI2; ANP32B
Species	Human
GeneID (Human)	10541
GeneID	10541
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ANP32B (Q92688)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92688
Protein Weight	29kDa
Shipping	Ice bag