


Anti-Anp32E Rabbit Polyclonal Antibody For Wb, Elisa - Q9Btt0

Item Number: CM0021346



Introduction

Anti-ANP32E rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat ANP32E. Targets the histone chaperone involved in H2A.Z removal. Suitable for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	ANP32E Rabbit pAb
Remarks/Alias	LANPL; LANP-L; ANP32E
Species	Human
Gene ID (Human)	81611
Gene ID	81611
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human ANP32E (NP_112182.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTT0
Protein Weight	31kDa
Shipping	Ice bag