

Anti-Anpep/Cd13 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P97449

Item Number: CM0011919

Introduction

Rabbit monoclonal antibody against ANPEP/CD13 (aminopeptidase N), validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studies in antigen presentation, peptide hormone processing, and amino acid transport.

[Learn More](#)

Parameter	Value
Product Name	ANPEP/CD13 Rabbit mAb
Remarks/Alias	Apn; AP-M; AP-N; Cd13; P150; CD13
Species	Mouse
Gene ID (Human)	290
Gene ID	16790
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 69-966 of mouse CD13. (NP_032512.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P97449
Protein Weight	109kDa
Shipping	Ice bag