

Anti-Anprb/Npr2 Rabbit Monoclonal Antibody For Wb, Elisa - P20594

Item Number: CM0002917

Introduction

Rabbit monoclonal antibody against human ANPRB/NPR2 (GC-B), validated for WB and ELISA. Cross-reacts with mouse. An essential tool for investigating cGMP signaling and skeletal growth pathways.

[Learn More](#)

Parameter	Value
Product Name	ANPRB/NPR2 Rabbit mAb
Remarks/Alias	GCB; AMD1; AMDM; ANPb; ECDM; GC-B; NPRB; SNSK; ANPRB; GUC2B; NPRBi; GUCY2B; ANPRB/NRP2
Species	Human
GeneID	4882
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 800-900 of human ANPRB/NRP2 (P20594).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P20594
Protein Weight	117kDa
Shipping	Ice bag