

Anti-Ant1/Ant2/Ant3/Ant4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P12235 / P05141 / P12236 / Q9H0C2

Item Number: CM0022252

Introduction

Rabbit polyclonal antibody targeting human ANT1/ANT2/ANT3/ANT4 (ADP/ATP translocases), suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial metabolism and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	ANT1/ANT2/ANT3/ANT4 Rabbit pAb
Remarks/Alias	ANT1/ANT2/ANT3/ANT4
Species	Human
GenelD (Human)	291/292/293/83447
GenelD	291/292/293/83447
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human ANT1 (NP_001142.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12235/P05141/P12236/Q9H0C2
Protein Weight	33kDa/32kDa/35kDa
Shipping	Ice bag