

Anti-Ant2/Slc25A5 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05141

Item Number: CM0006737

Introduction

Recombinant rabbit monoclonal antibody against human ANT2/SLC25A5 (ADP/ATP translocase 2), validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for mitochondrial research and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ANT2/SLC25A5 Rabbit mAb
Remarks/Alias	T2; T3; 2F1; AAC2; ANT2; ANT2/SLC25A5
Species	Human
GeneID (Human)	292
GeneID	292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-220 of human ANT2/SLC25A5 (NP_001143.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05141
Protein Weight	33kDa
Shipping	Ice bag