


Anti-Anxa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04083

Item Number: CM0017211



Introduction

Rabbit polyclonal antibody targeting human ANXA1 (Annexin A1), suitable for WB, IHC-P, IF/ICC, and ELISA. Detects human and rat ANXA1; ideal for inflammation and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	ANXA1 Rabbit pAb
Remarks/Alias	ANX1; LPC1; ANXA1
Species	Human
GeneID (Human)	301
GeneID	301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human ANXA1 (NP_000691.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P04083
Protein Weight	39kDa
Shipping	Ice bag