


Anti-Anxa13 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P27216

Item Number: CM0014752



Introduction

Rabbit polyclonal antibody targeting human annexin A13 (ANXA13), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat ANXA13. Suitable for studying calcium-dependent membrane binding and intestinal epithelial biology.

[Learn More](#)

Parameter	Value
Product Name	ANXA13 Rabbit pAb
Remarks/Alias	ISA; ANX13; ANXA13
Species	Human
Gene ID (Human)	312
Gene ID	312
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-316 of human ANXA13 (NP_004297.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27216
Protein Weight	35kDa
Shipping	Ice bag