

Anti-Anxa3 Polyclonal Antibody For Wb And Elisa - P12429

Item Number: CM0017409



Introduction

Rabbit polyclonal antibody targeting human Annexin A3 (ANXA3) validated for Western blot and ELISA, with reactivity in human and rat. Ideal for studying phospholipase A2 inhibition and coagulation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ANXA3 Rabbit pAb
Remarks/Alias	ANX3, ANXA3
Species	Human
GeneID	306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human ANXA3 (NP_005130.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P12429
Protein Weight	36kDa
Shipping	Ice bag