

Anti-Anxa3 Polyclonal Antibody For Wb - P12429

Item Number: CM0025389



Introduction

Rabbit polyclonal antibody against human ANXA3, validated for WB and ELISA. Cross-reacts with mouse and rat. Targets the phospholipase A2 inhibitor annexin A3, suitable for coagulation and phospholipid signaling studies.

[Learn More](#)

Parameter	Value
Product Name	ANXA3 Rabbit pAb
Remarks/Alias	ANX3; ANXA3
Species	Human
GeneID (Human)	306
GeneID	306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human ANXA3 (NP_005130.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12429
Protein Weight	36kDa
Shipping	Ice bag