

Anti-Anxa7 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P20073

Item Number: CM0003812



Introduction

High-quality ANXA7 Rabbit Monoclonal Antibody for WB, IF/ICC, ELISA. Detects human, mouse, rat ANXA7. Ideal for studying annexin A7's role in calcium-mediated exocytosis and membrane fusion.

[Learn More](#)

Parameter	Value
Product Name	ANXA7 Rabbit mAb
Remarks/Alias	SNX; ANX7; SYNEXIN; ANXA7
Species	Human
Gene ID (Human)	310
Gene ID	310
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ANXA7 (NP_001147.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20073
Protein Weight	53kDa
Shipping	Ice bag