

Anti-Anxa8L2 Rabbit Pab For Wb And Elisa - Q5Vt79

Item Number: CM0026169

Introduction

Rabbit polyclonal antibody targeting annexin A8-like protein 1 (ANXA8L2, Q5VT79). Validated for WB and ELISA, reactive with human, mouse, and rat samples. Suitable as IVD raw material for diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ANXA8L2 Rabbit pAb
Remarks/Alias	ANXA8; ANXA8L2; VAC-beta; bA145E20.2
Species	Human
Gene ID (Human)	728113
Gene ID	728113
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-276 of human ANXA8L2 (NP_001265853.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VT79
Protein Weight	37kDa
Shipping	Ice bag