


Anti-Aoc3 Rabbit Polyclonal Antibody For Wb And Elisa - Q16853

Item Number: CM0023158



Introduction

Rabbit polyclonal antibody against human amine oxidase copper-containing 3 (AOC3/VAP-1). Suitable for Western blotting and ELISA; cross-reacts with mouse and rat. Ideal for studies of amine oxidase function and lymphocyte adhesion.

[Learn More](#)

Parameter	Value
Product Name	AOC3 Rabbit pAb
Remarks/Alias	HPAO; SSAO; VAP1; VAP-1; AOC3
Species	Human
Gene ID (Human)	8639
Gene ID	8639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-190 of human AOC3 (NP_003725.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16853
Protein Weight	85kDa
Shipping	Ice bag