


Anti-Aox1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q06278

Item Number: CM0027923



Introduction

High-quality rabbit polyclonal antibody against human AOX1 (aldehyde oxidase 1). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for xenobiotic metabolism and drug metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AOX1 Rabbit pAb
Remarks/Alias	AO; AOH1; AOX1
Species	Human
GeneID (Human)	316
GeneID	316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-590 of human AOX1 (NP_001150.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06278
Protein Weight	148kDa
Shipping	Ice bag