


Anti-Ap1B1 Rabbit Polyclonal Antibody For Wb, Elisa - Q10567

Item Number: CM0020308



Introduction

Rabbit polyclonal antibody against AP1B1, validated for Western blotting and ELISA. Reacts with human, mouse, and rat AP1B1, a clathrin adaptor subunit involved in protein sorting at the TGN and endosomes.

[Learn More](#)

Parameter	Value
Product Name	AP1B1 Rabbit pAb
Remarks/Alias	ADTB1; BAM22; KIDAR; AP105A; CLAPB2; AP1B1
Species	Human
GeneID (Human)	162
GeneID	162
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human AP1B1 (NP_001118.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q10567
Protein Weight	105kDa
Shipping	Ice bag