

Anti-AP1G1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O43747

Item Number: CM0003510



Introduction

Anti-AP1G1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes AP-1 complex subunit gamma-1 in human, mouse, and rat samples. Suitable for studying clathrin-mediated protein sorting and trafficking.

[Learn More](#)

Parameter	Value
Product Name	AP1G1 Rabbit mAb
Remarks/Alias	ADTG; CLAPG1; USRISD; AP1G1
Species	Human
Gene ID (Human)	164
Gene ID	164
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 723-822 of human AP1G1 (O43747)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43747
Protein Weight	91kDa
Shipping	Ice bag