

Anti-Ap1M1 Polyclonal Antibody For Wb And Elisa - Q9Bxs5

Item Number: CM0015609

Introduction

AP1M1 Rabbit pAb, a polyclonal antibody for WB and ELISA, detects AP-1 complex subunit mu-1 in human, mouse, and rat. Used in protein sorting and clathrin-mediated trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	AP1M1 Rabbit pAb
Remarks/Alias	AP47, mu1A, CLTNM, MU-1A, CLAPM2, AP1M1
Species	Human
Gene ID (Human)	8907
Gene ID	8907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human AP1M1 (NP_115882.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXS5
Protein Weight	49kDa
Shipping	Ice bag