

Anti-Ap1M2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6Q5

Item Number: CM0029657

Introduction

Rabbit polyclonal antibody targeting AP1M2 (mu-2 subunit of AP-1 complex). Suitable for WB, IF/ICC, ELISA; reacts with human, mouse, rat. For research in protein sorting and clathrin-mediated trafficking.

[Learn More](#)

Parameter	Value
Product Name	AP1M2 Rabbit pAb
Remarks/Alias	mu2; MU1B; MU-1B; HSMU1B; AP1-mu2; AP1M2
Species	Human
GeneID (Human)	10053
GeneID	10053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 164-423 of human AP1M2 (NP_005489.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6Q5
Protein Weight	48kDa
Shipping	Ice bag