

Anti-Ap2A1 Polyclonal Antibody For Wb, Elisa - O95782

Item Number: CM0025264

Introduction

High-quality rabbit polyclonal antibody targeting human AP2A1 (adaptor protein complex 2 subunit alpha-1). Validated for Western blot and ELISA, cross-reactive with mouse and rat. Ideal for endocytosis and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	AP2A1 Rabbit pAb
Remarks/Alias	ADTAA; CLAPA1; AP2-ALPHA; AP2A1
Species	Human
GeneID (Human)	160
GeneID	160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AP2A1 (NP_055018.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95782
Protein Weight	108kDa
Shipping	Ice bag