

Anti-Ap2A2/Alpha2-Adaptin Monoclonal Antibody For Wb, Elisa - O94973

Item Number: CM0003898

Introduction

CamelBio's Alpha2-adaptin/AP2A2 rabbit monoclonal antibody is validated for WB and ELISA. Detects human and mouse AP2A2, a key clathrin adaptor in endocytosis. Ideal for neuroscience and cell biology research.

[Learn More](#)

Parameter	Specification
Product Name	Alpha2-adaptin/AP2A2 Rabbit mAb
Remarks/Alias	HIP9; HYPJ; ADTAB; HIP-9; CLAPA2; Alpha2-adaptin/AP2A2
Species	Human
Gene ID (Human)	161
Gene ID	161
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 590-939 of human Alpha2-adaptin/AP2A2 (NP_036437.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94973
Protein Weight	104kDa
Shipping	Ice bag