

Anti-AP2M1 Monoclonal Antibody For Wb And Elisa - Q96Cw1

Item Number: CM0002165



Introduction

Rabbit monoclonal antibody specific to human AP2M1 (AP-2 mu subunit), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Detects a 50kDa protein. Ideal for endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	AP2M1 Rabbit mAb
Remarks/Alias	mu2; AP50; MRD60; CLAPM1; AP2M1
Species	Human
Gene ID (human)	1173
Gene ID	1173
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 300-400 of human AP2M1 (Q96CW1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96CW1
Protein Weight	50kDa
Shipping	Ice bag