

Anti-Ap2M1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Cw1

Item Number: CM0016197



Introduction

Rabbit polyclonal antibody targeting AP2M1, a key subunit of the AP-2 adaptor complex involved in clathrin-mediated endocytosis. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	AP2M1 Rabbit pAb
Remarks/Alias	mu2; AP50; MRD60; CLAPM1; AP2M1
Species	Human
GeneID (Human)	1173
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-433 of human AP2M1 (NP_004059.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96CW1
Protein Weight	50kDa
Shipping	Ice bag