

Anti-AP2M1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Cw1

Item Number: CM0019136



Introduction

CamelBio's AP2M1 rabbit polyclonal antibody, validated for WB and ELISA, targets the AP-2 complex mu subunit involved in clathrin-mediated endocytosis. Suitable for human samples, with a detected molecular weight of 50kDa.

[Learn More](#)

Parameter	Value
Product Name	AP2M1 Rabbit pAb
Remarks/Alias	mu2; AP50; MRD60; CLAPM1; AP2M1
Species	Human
GeneID (Human)	1173
GeneID	1173
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96CW1
Protein Weight	50kDa
Shipping	Ice bag