


Anti-Ap2S1 Rabbit Monoclonal Antibody For Wb And Elisa - P53680

Item Number: CM0010806



Introduction

Anti-AP2S1 rabbit monoclonal antibody, 17 kDa, validated for Western blot (WB) and ELISA. Cross-reacts with mouse; ideal for studying clathrin-mediated endocytosis and adaptor protein complex 2 function.

[Learn More](#)

Parameter	Value
Product Name	AP2S1 Rabbit mAb
Remarks/Alias	AP17; FBH3; HHC3; FBHOk; CLAP52; AP2S1
Species	Human
GeneID	1175
Immunogen	Synthetic Peptide, A synthetic peptide corresponding to a sequence within amino acids 1-100 of human AP2S1 (P53680)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P53680
Protein Weight	17kDa
Shipping	Ice bag