

Anti-Ap3D1 Rabbit Polyclonal Antibody For Wb, Elisa - O14617

Item Number: CM0018013



Introduction

Anti-AP3D1 Rabbit Polyclonal Antibody validated for WB and ELISA. Cross-reacts with human and mouse AP3D1. Ideal for research on AP-3 adaptor complex, lysosomal trafficking, and Hermansky-Pudlak syndrome.

[Learn More](#)

Parameter	Value
Product Name	AP3D1 Rabbit pAb
Remarks/Alias	ADTD; HPS10; hBLVR; AP3D1
Species	Human
GeneID (Human)	8943
GeneID	8943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-700 of human AP3D1 (NP_003929.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14617
Protein Weight	130kDa
Shipping	Ice bag