


Anti-Apaf1 Rabbit Monoclonal Antibody For Wb And Elisa - O14727

Item Number: CM0008325



Introduction

Recombinant rabbit monoclonal antibody targeting human APAF1. Validated for WB and ELISA, cross-reacts with human and mouse. Detects endogenous APAF1 protein (38-142 kDa).

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] APAF1 Rabbit mAb
Remarks/Alias	CED4, APAF-1
Species	Human
Gene ID (Human)	317
Gene ID	317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-600 of human APAF1 (NP_863651.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14727
Protein Weight	38-142 kDa
Shipping	Ice bag