

Anti-Apba1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q02410

Item Number: CM0023631



Introduction

Rabbit polyclonal antibody targeting APBA1 (MINT1/X11). Validated for IHC-P, IF/ICC, ELISA in human and mouse tissues. Essential tool for studying synaptic vesicle exocytosis, APP processing, and NMDA receptor trafficking. Reliable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	APBA1 Rabbit pAb
Remarks/Alias	X11; X11A; LIN10; MINT1; D9S411E; X11ALPHA; APBA1
Species	Human
Gene ID (Human)	320
Gene ID	320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-390 of human APBA1 (NP_001154.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q02410
Protein Weight	93kDa
Shipping	Ice bag