

Anti-Apbb1 Rabbit Polyclonal Antibody For Wb, Elisa - 000213

Item Number: CM0024557



Introduction

Rabbit polyclonal antibody targeting human APBB1 (FE65) for WB and ELISA applications. Cross-reacts with mouse. Ideal for studying DNA damage response, transcription regulation, and Alzheimer's disease pathways.

[Learn More](#)

Parameter	Value
Product Name	APBB1 Rabbit pAb
Remarks/Alias	RIR; FE65; MGC:9072; APBB1
Species	Human
Gene ID (Human)	322
Gene ID	322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 359-708 of human APBB1 (NP_663722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00213
Protein Weight	77kDa
Shipping	Ice bag