


## Anti-Apc Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P25054

Item Number: CM0025687



### Introduction

Anti-APC polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat APC (312 kDa). Suitable for studying tumor suppression, Wnt signaling, and cell migration.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | APC Rabbit pAb                                                                                                                           |
| Remarks/Alias    | GS; DP2; DP3; BTPS2; DESMD; DP2.5; PPP1R46; APC                                                                                          |
| Species          | Human                                                                                                                                    |
| GeneID (human)   | 324                                                                                                                                      |
| GeneID           | 324                                                                                                                                      |
| Immunogen        | Recombinant protein. Recombinant fusion protein containing a sequence corresponding to amino acids 2600-2760 of human APC (NP_000029.2). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                                    |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                        |
| SWISS            | P25054                                                                                                                                   |
| Protein Weight   | 312kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                  |