


Anti-Apc Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P25054

Item Number: CM0005625



Introduction

Rabbit monoclonal antibody against human APC (adenomatous polyposis coli) protein, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for tumor suppressor and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit mAb
Remarks/Alias	GS; DP2; DP3; BTPS2; DESMD; DP2.5; PPP1R46; APC
Species	Human
GeneID	324
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2600-2760 of human APC (NP_000029.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25054
Protein Weight	312kDa
Shipping	Ice bag