

Anti-Apc Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P25054

Item Number: CM0002567



Introduction

Rabbit monoclonal antibody against human APC (adenomatous polyposis coli), a tumor suppressor in Wnt signaling. Validated for WB, IP, ELISA; cross-reacts with mouse and rat. For cancer research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit mAb
Remarks/Alias	GS; DP2; DP3; BTP52; DESMD; DP2.5; PPP1R46; APC
Species	Human
GeneID (Human)	324
GeneID	324
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human APC (NP_000029.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25054
Protein Weight	312kDa
Shipping	Ice bag