


Anti-Apcs Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02743

Item Number: CM0023156



Introduction

High-quality rabbit polyclonal antibody targeting human APCS (Serum amyloid P-component). Validated for WB, IHC-P, ELISA in human, mouse, rat samples. Suitable for diagnostic assay development and research.

[Learn More](#)

Parameter	Value
Product Name	APCS Rabbit pAb
Remarks/Alias	SAP; PTX2; HEL-S-92n; APCS
Species	Human
Gene ID (Human)	325
Gene ID	325
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-223 of human APCS (NP_001630.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02743
Protein Weight	25kDa
Shipping	Ice bag