

Anti-Apeh Polyclonal Antibody For Wb, If/Icc, Elisa - P13798

Item Number: CM0018435

Introduction

Rabbit polyclonal antibody targeting human APEH (Acylamino-acid-releasing enzyme), suitable for WB, IF/ICC, ELISA applications. Cross-reacts with mouse and rat. Useful for studying N-acylated peptide degradation and oxidized protein hydrolase activity.

[Learn More](#)

Parameter	Value
Product Name	APEH Rabbit pAb
Remarks/Alias	APH; OPH; AARE; ACPH; D3S48E; D3F15S2; DNF15S2; APEH
Species	Human
GeneID (Human)	327
GeneID	327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human APEH (NP_001631.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13798
Protein Weight	81kDa
Shipping	Ice bag