

Anti-Apeh Polyclonal Antibody For Wb, Ip, Elisa - P13798

Item Number: CM0025042

Introduction

Anti-APEH Rabbit Polyclonal Antibody validated for WB, IP, and ELISA applications, with cross-reactivity in human, mouse, and rat. Targets acylamino-acid-releasing enzyme (P13798), an 81 kDa cytoplasmic hydrolase for N-acylated peptide degradation. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] APEH Rabbit pAb
Remarks/Alias	APH; OPH; AARE; ACPH; D3S48E; D3F15S2; DNF15S2; EH
Species	Human
Gene ID (Human)	327
Gene ID	327
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human APEH (NP_001631.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13798
Protein Weight	81kDa
Shipping	Ice bag