


Anti-Apex1/Ape1 Polyclonal Antibody For Wb, If/Icc, Elisa - P27695

Item Number: CM0015642



Introduction

Rabbit polyclonal antibody targeting human APEX1/APE1 (SWISS P27695), suitable for WB, IF/ICC, and ELISA. Recognizes the 36kDa DNA repair nuclease/redox regulator, ideal for BER pathway and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	APEX1/APE1 Rabbit pAb
Remarks/Alias	APE; APX; APE1; APEN; APEX; HAP1; REF1; APEX1/APE1
Species	Human
GeneID (Human)	328
GeneID	328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human APEX1/APE1 (NP_542380.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P27695
Protein Weight	36kDa
Shipping	Ice bag