


Anti-Apex1/Ape1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Ip, Elisa - P27695

Item Number: CM0026130



Introduction

KO validated rabbit polyclonal antibody to human APEX1/APE1 (AP endonuclease 1, REF1). Detects a 36 kDa band. Ideal for WB, IF/ICC, IP, ELISA in oxidative stress and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] APEX1/APE1 Rabbit pAb
Remarks/Alias	APE; APX; APE1; APEN; APEX; HAP1; REF1; E1
Species	Human
Gene ID (Human)	328
Gene ID	328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human APEX1/APE1 (NP_542380.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P27695
Protein Weight	36kDa
Shipping	Ice bag