


Anti-Aph1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Bi3

Item Number: CM0018321



Introduction

APH1A Rabbit pAb is a rabbit polyclonal antibody targeting the gamma-secretase subunit APH1A. Suitable for WB, IF/ICC, ELISA with human and mouse reactivity. Supports Notch and APP cleavage pathway research.

[Learn More](#)

Parameter	Value
Product Name	APH1A Rabbit pAb
Remarks/Alias	APH-1; APH-1A; CGI-78; 6530402N02Rik; APH1A
Species	Human
GeneID (Human)	51107
GeneID	51107
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human APH1A (NP_001071096.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96Bi3
Protein Weight	29kDa
Shipping	Ice bag