

Anti-Apip Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Gx9

Item Number: CM0028244



Introduction

High-quality rabbit polyclonal antibody against human APIP (methylthioribulose-1-phosphate dehydratase). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for studying methionine salvage and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	APIP Rabbit pAb
Remarks/Alias	APIP2; CGI29; hAPIP; CGI-29; MMRP19; APIP
Species	Human
GeneID (Human)	51074
GeneID	51074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human APIP (NP_057041.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GX9
Protein Weight	27kDa
Shipping	Ice bag