

Anti-Apitd1 Rabbit Pab - Q8N2Z9

Item Number: CM0025855



Introduction

Polyclonal antibody against human APITD1 (CENPS), a key player in Fanconi anemia pathway and kinetochore assembly. Validated for WB and ELISA.

[Learn More](#)

Parameter	Value
Product Name	APITD1 Rabbit pAb
Remarks/Alias	MHF1; APITD1; CENP-S; FAAP16
Species	Human
GeneID (Human)	378708
GeneID	378708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human APITD1 (NP_954988.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N2Z9
Protein Weight	16kDa
Shipping	Ice bag