

Anti-Aplf Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iw19

Item Number: CM0020106



Introduction

Anti-APLF rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects human APLF with rat cross-reactivity. Ideal for DNA damage repair, histone chaperone, and NHEJ research. 57 kDa.

[Learn More](#)

Parameter	Value
Product Name	APLF Rabbit pAb
Remarks/Alias	APFL; PALF; Xip1; ZCCHH1; C2orf13; APLF
Species	Human
GeneID (Human)	200558
GeneID	200558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human APLF (NP_775816.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8IW19
Protein Weight	57kDa
Shipping	Ice bag