


Anti-Aplp1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P51693

Item Number: CM0025933



Introduction

Anti-APLP1 rabbit polyclonal antibody validated for WB, IF-P, and ELISA applications. Detects APLP1 in human, mouse, and rat samples. Ideal for studying postsynaptic function, neuronal apoptosis, and Alzheimer's-related pathways.

[Learn More](#)

Parameter	Value
Product Name	APLP1 Rabbit pAb
Remarks/Alias	APLP; APLP1
Species	Human
GeneID (Human)	333
GeneID	333
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-580 of human APLP1 (NP_005157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51693
Protein Weight	72kDa
Shipping	Ice bag