


Anti-Aplp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51693

Item Number: CM0026913



Introduction

Rabbit polyclonal antibody targeting human APLP1 (Amyloid beta precursor like protein 1). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for neuroscience and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	APLP1 Rabbit pAb
Remarks/Alias	APLP; APLP1
Species	Human
Gene ID (Human)	333
Gene ID	333
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 285-490 of human APLP1 (NP_005157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51693
Protein Weight	72kDa
Shipping	Ice bag