

Anti-Aplp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q06481

Item Number: CM0030576



Introduction

High-quality APLP2 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA, cross-reacting with human, mouse, and rat. Ideal for amyloid and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	APLP2 Rabbit pAb
Remarks/Alias	APPH, APPL2, CDEBP, APLP-2, APLP2
Species	Human
Gene ID (Human)	334
Gene ID	334
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-580 of human APLP2 (NP_001633.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06481
Protein Weight	87kDa
Shipping	Ice bag