


Anti-Apoa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02647

Item Number: CM0014824



Introduction

Anti-APOA1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human APOA1 with cross-reactivity to mouse and rat. Ideal for lipid metabolism, atherosclerosis, and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	APOA1 Rabbit pAb
Remarks/Alias	HPALP2; apo(a); APOA1
Species	Human
Gene ID (Human)	335
Immunogen	Recombinant fusion protein containing amino acids 25-267 of human APOA1 (NP_000030.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P02647
Protein Weight	28kDa
Shipping	Ice bag