

Anti-Apoa4 Polyclonal Antibody For Wb, Elisa - P06727

Item Number: CM0028556



Introduction

Rabbit polyclonal antibody targeting human APOA4, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for lipid metabolism and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	APOA4 Rabbit pAb
Remarks/Alias	APOA4
Species	Human
GeneID (Human)	337
GeneID	337
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 227-396 of human APOA4 (NP_000473.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06727
Protein Weight	45kDa
Shipping	Ice bag