

Anti-Apoa4 Rabbit Monoclonal Antibody For Wb, Elisa - P06727

Item Number: CM0007370



Introduction

Recombinant rabbit monoclonal antibody against human APOA4. Validated for WB and ELISA. Ideal for lipid metabolism research, cholesterol transport studies, and diagnostic development. Works on human samples. High-affinity antibody targeting amino acids 21-396. Recognizes a ~45kDa protein.

[Learn More](#)

Parameter	Value
Product Name	APOA4 Rabbit mAb
Remarks/Alias	APOA4
Species	Human
Gene ID (Human)	337
Gene ID	337
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 21-396 of human APOA4 (NP_000473.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06727
Protein Weight	45kDa
Shipping	Ice bag