

## Anti-Apoa5 Polyclonal Antibody For Wb, Elisa - Q6Q788

Item Number: CM0018252

### Introduction

Rabbit polyclonal antibody targeting human APOA5, validated for Western blotting and ELISA. Recognizes human and mouse APOA5, a key regulator of triglyceride metabolism. Suitable for lipid research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | APOA5 Rabbit pAb                                                                                                   |
| Remarks/Alias    | RAP3; APOAV; APOA5                                                                                                 |
| Species          | Human                                                                                                              |
| Gene ID          | 116519 (Human)                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 207-366 of human APOA5 (NP_443200.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q6Q788                                                                                                             |
| Protein Weight   | 41kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |